


## (R)EVOLUTION®

The Champions Implant (R)Evolution® is a screw-shaped two-piece implant that is inserted to bone level. The implant has diameters of 3.5/4.0/4.5/5.5 mm. It is available in lengths of 6.5/8/10/12/14/16 mm. The implant covers a wide range of indications. The implant is suitable for the insertion in the long-term healed jaw bone (delayed implantation), also for the delayed insertion (6-8 weeks after the tooth extraction) and under certain conditions for immediate implantation (implantation directly after the tooth extraction). The Champions Implant (R)Evolution® can be placed in bone in the maxilla and mandible (D1-D4) with and without augmentations.

Please note the restrictions of indications in the Instructions for Use.

Material: Titanium grade 4



Diameter	Length	Height	Art. No.	Type
3.5 mm	6.5 mm	3.5 mm	0800	●
	8 mm	3.5 mm	0815	
	10 mm	3.5 mm	0830	
	12 mm	3.5 mm	0845	
	14 mm	3.5 mm	0860	
	16 mm	3.5 mm	0875	
4.0 mm	6.5 mm	3.5 mm	0802	
	8 mm	3.5 mm	0817	
	10 mm	3.5 mm	0832	
	12 mm	3.5 mm	0847	
	14 mm	3.5 mm	0862	
	16 mm	3.5 mm	0877	
4,5 mm	6.5 mm	3.5 mm	0805	
	8 mm	3.5 mm	0820	
	10 mm	3.5 mm	0835	
	12 mm	3.5 mm	0850	
	14 mm	3.5 mm	0865	
	16 mm	3.5 mm	0880	
5.5 mm	6.5 mm	3.5 mm	0810	
	8 mm	3.5 mm	0825	
	10 mm	3.5 mm	0840	
	12 mm	3.5 mm	0855	
	14 mm	3.5 mm	0870	
	16 mm	3.5 mm	0885	